

Table 1: Primers and Probes used for One step RT-PCR of vesicular disease viruses

		Oligo name	Sequence (5' – 3')	Working conc (pmols/ul)
FMDV ^{5UTR}	Forward Primer	SA-IR-219-246F	CAC YTY AAG RTG ACA YTG RTA CTG GTA C	10
	Reverse Primer	SA-IR-315-293R	CAG ATY CCR AGT GWC ICI TGT TA	10
	Taqman PROBE	SAmulti2-P-IR-292-269R	CCT CGG GGT ACC TGA AGG GCA TCC	5
FMDV ^{3D}	Forward Primer	3DF	ACT GGG TTT TAC AAA CCT GTG A	10
	Reverse Primer	3DR	GCG AGT CCT GCC ACG GA	10
	Taqman PROBE	3DP	TCC TTT GCA CGC CGT GGG AC	5